



beurer
health and well-being

Clear airways for the littlest ones: Beurer expands BabyCare range

Ulm, 17 November 2019 – Coughing, sniffles, hoarseness: if your child is ill, it can be a very stressful time for both the children and parents. This is where products which not only provide fast relief but are also easy to use come into their own. At MEDICA, Beurer is presenting a nasal aspirator and two nebulisers in the Kids version, providing clear airways for babies and small children.

NA 20 – nasal aspirator: effective, gentle and easy to use

The **NA 20** nasal aspirator allows for gentle cleaning of the nostrils for babies or small children, making it easier for them to breathe, sleep and eat. It is ideal for use during the colder months. The product is suitable for removing liquid to solid nasal secretions. The ergonomic design, handy one-button operation and controlled suction pressure (maximum 66 kPa) make the device especially easy to use. The NA 20 is quiet and battery operated. The removable elements are easy to clean. Two soft silicone attachments which can be disinfected round out the items included in delivery.

Effective inhalation from the 1st month of life: IH 26 and IH 58 Kids

The **IH 26** and the **IH 58** in the Kids edition are inhalation devices for atomising liquid medication using compressed-air technology. The medical devices are suitable for treating the upper and lower airways. In doing so, the application relieves the side effects of colds and speeds up recovery. However, it can also be used for other respiratory diseases such as asthma. Both nebulisers are characterised by their child-friendly design. The attachable jolly giraffe will help distract the child in a playful way during treatment. The high nebulisation performance of approx. 0.3 millilitres per minute for the IH 26 Kids, and 0.25 millilitres per minute for the IH 58 Kids, allows inhalation to be carried out effectively in next to no time. The devices can be disinfected, are easy to operate and come with the following accessories: silicone children's and baby mask, angled fitting, mouthpiece, atomiser, 10 filters, hose.

Beurer GmbH
Germany, **Ulm**
Soeflinger Strasse 218
89077 Ulm
Germany

Contact
Anna-Katharina Peuker
presse@beurer.de

Website
www.beurer.com
www.beurer-healthguide.com

Social Media:



Technical data:



IH 26 Kids nebuliser

- Compressed-air nebuliser
- Particle size (MMAD) 3.07 μm
- Working pressure approx. 0.8-1.45 bar
- Nebulisation performance approx. 0.3 ml/min
- Incl. extensive child-friendly accessories



IH 58 Kids nebuliser

- Compressed-air nebuliser
- Particle size (MMAD) 4.12 μm
- Working pressure approx. 0.25-0.5 bar
- Nebulisation performance approx. 0.25 ml/min
- Quiet and especially practical when travelling
- Incl. extensive child-friendly accessories

All of the highlights will be on show at the Beurer MEDICA stand from **18 to 21 November 2019** in **Hall 10, Booth D20**.

About Beurer

Beurer was founded in 1919 in Ulm and is synonymous with health and well-being. Today the long-standing company leads the way in several product fields in this segment; the company is the market leader in Europe in the area of flexible heating and the market leader in Germany for blood pressure monitors and massage products. Beurer is also one of the leading European suppliers of personal scales. This portfolio is constantly being developed further in all areas and offers products for use at home. It includes personal scales, kitchen scales, luggage scales, air cleaners, air dehumidifiers, air washers, thermo hygrometers, aroma diffusers, blood pressure monitors, blood glucose monitors, mobile ECG devices, hearing amplifiers, pulse oximeters, nebulisers, clinical thermometers, a sleep sensor, wake-up lights, a snore stopper, brightlights, infrared lamps, TENS/EMS devices, massage products (foot massage, Shiatsu massage, massage chairs, massage cushions, fasciae massage, compression leg therapy), a relaxing aid, hair removal devices (classic, IPL, laser), FaceCare and HairCare products, a cellulite releaZer, cosmetics mirrors, manicure/pedicure sets, a BabyCare line, basal thermometers, activity sensors, an activity watch and heart rate monitors. USB ports and *Bluetooth*[®] enable an increasing number of Connect products to connect to the growing Beurer app world and/or software. The family-run company operates a global distribution network in more than 100 countries and currently has a workforce of around 900. Further information is available at www.beurer.com.