

Press release

## Optimum monitoring of blood glucose levels thanks to *Bluetooth*<sup>®</sup> function

Ulm, 20 November 2024 – The GL 60 Blood Glucose Monitor from Beurer sets the benchmark when it comes to blood glucose monitoring and provides a reliable and intuitive way for users to check their blood glucose level. It can also be used by people who have gestational diabetes. Boasting a wealth of innovative features, the *Bluetooth*<sup>®</sup>-enabled GL 60 is easy to use and plays a valuable role in the provision of state-of-the-art healthcare in line with each person's individual needs.

### Maximum precision and safety for the user

The device requires a minimum of 0.5 µl of blood to ensure it can be used safely and precisely. Before the measurement is taken, the device automatically checks that the sample is sufficient. Safety is further boosted by the wide haematocrit range from 0 to 70%, which also makes the device suitable for monitoring gestational diabetes. The wide test strips ensure the device is easy to use. The GL 60 has been clinically tested and is certified by the Institute for Diabetes Karlsburg.

### Advanced warning functions and mobile integration into the “beurer HealthManager Pro” app

The GL 60 also features an integrated ketone and glucose warning to inform users directly about potential health risks. The measured values can be transferred directly to the “beurer HealthManager Pro” app via *Bluetooth*<sup>®</sup> or USB and evaluated there. Moreover, the device is compatible with the app as well as PC software. This enables health data to be monitored in real time and exported, so users can easily share it with doctors or keep it for their own records.


### Comprehensive health management with “beurer HealthManager Pro”

Alongside blood glucose monitoring, the “beurer HealthManager Pro” app offers comprehensive digital health management. The app brings together numerous health parameters such as blood glucose, blood pressure, weight, oxygen saturation, activity and temperature. As a certified medical device, it meets the high requirements of the relevant European directives and provides a detailed overview of users' health data.



## Advantages of the “beurer HealthManager Pro” app at a glance:

- **Medical device:**  
Complies with European directives and provides detailed measurement values.
- **Always available:**  
Measurements and progress graphs can be viewed anywhere and at any time.
- **User-friendly:**  
Modern design and intuitive operation.
- **Accessible:** Usable by all and easy to understand.
- **Export function:** Quick and easy transfer to doctors.



The GL 60 allows you to **calculate the average value separately** for periods of 7, 14, 30 and 90 days. This function enables all measured values, or values that have been marked specifically, to be evaluated so that users can obtain a comprehensive, precise analysis. Optimum diabetes management is ensured by the measurement marking: before eating, after eating, general marking.

The GL 60 is battery-powered and supports **several languages** (German, English, Spanish, French, Italian and Turkish). An illuminated insertion slot makes it effortless to use, even in poor lighting conditions. The GL 60's wide test strips and easy-to-read LCD display with large digits make it the ideal choice for everyday use.

### Beurer GmbH

Söflinger Straße 218  
D- 89077 Ulm  
Tel. +49 731/3989-0  
Fax +49 731/3989-295  
[www.beurer.com](http://www.beurer.com)

### Beurer Press Contact

Simon Schwörer  
E-Mail HQ: [presse@beurer.de](mailto:presse@beurer.de)

### Social Media



## **About Beurer**

Beurer was founded in 1919 in Ulm and is synonymous with health and well-being. Today, the family-run company leads the way in several areas of this segment. Originally a manufacturer of heat pads and heated underblankets, the company has grown into a full-range supplier offering more than 500 products. Beurer offers various products for use at home in the categories of health, wellbeing, beauty, personal care, home, fitness and baby, including blood pressure and blood glucose monitors, clinical thermometers, pulse oximeters, nebulisers, hearing amplifiers, daylight therapy lamps, insect bite healers, a mobile ECG device and an emergency wristband, as well as products for pain therapy and women's health. The range also covers the areas of flexible heating, weight & diagnosis, air & aroma, massage & shiatsu and sleep & rest. What's more, Beurer offers various beauty products in the categories of FaceCare, BodyCare, HairCare, hair removal, manicure & pedicure, as well as a men's grooming range and electric toothbrushes in the personal care category. The company offers a range of fitness products in its Antelope by beurer and PhysioLine product groups, along with numerous abdominal muscle trainers. The baby category contains special products adapted to the needs of children, babies and their parents. With its "Connect" product group, the health specialist offers a wide range of apps that can be combined with Beurer products and make everyday life easier. Around 1,700 employees worldwide and a global distribution network spanning more than 100 countries ensure the company's ever-growing success. You can find more information at [www.beurer.com](http://www.beurer.com).