



OT 30
BasalControl
Smart

beurer
Woman's Life

BasalControl Smart

Ovulation thermometer



Can be used with the app "Ovy" - Data transmission of temperature with *Bluetooth*[®] low energy technology



Calculation of the cycle for natural family planning



Documentation of more body signals



- Digital ovulation checking thermometer
- Measurement of basal body temperature
- Display of 2 digits after comma
- For oral and vaginal measurement
- Illuminated LED ring signals *Bluetooth*[®] transmission
- Measurement accuracy:
± 0.05 °C: 35.00 °C - 38.00 °C at ambient temperature from 18 °C to 28 °C;
outside this range ±0.1 °C
- Measurement range: 32.00 - 43.99 °C
- Flexible measurement tip
- Waterproof
- Can be disinfected
- 30 memory spaces
- Signal tone at end of measurement (On/Off)
- With protection cap
- Including 1x Battery CR2032
- Battery changing indicator
- Automatic switch off
- Indication of °C/°F
- Product measurements (LxWxH): 105 x 33 x 17.6 mm
- Product weight: approx. 20.1 g (without battery)
- Medical device
- 5 year guarantee
- Sales unit: 4 / Shipping carton: 12 x 4
- EAN no.: 4211125 79106 3
- Item no.: 791.06



OT 30
BasalControl
Smart

beurer
Woman's Life

- Pregnancy planning or cycle tracking
- Record temperature, body signals and interference factors in the app
- Ovy Care provides tips on nutrition, fitness and beauty according to the current cycle phase
- Medical device class I



Advanced cycle monitoring in just four steps


Wake up and take basal body temperature

1




2

Transfer temperature to the app



Record body signals in the app

3



4

Prediction of fertile and infertile days



List of system requirements and compatible devices

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.