

## BC 28 Wrist Blood Pressure Monitor

- Risk indicator: Classification of measurements with coloured scale
- Arrhythmia detection: Warning in case of possible heart rhythm disturbances
- Ø morning and evening blood pressure for the last 7 days
- Two user memories for 60 measured values each





## **BM 28** Upper Arm Blood Pressure Monitor

- Resting indicator\* for exact measurement results
- Cuff position control: For ensuring the cuff is correctly positioned on the upper arm



\*Patent EP 1 673 009



## **BM 54 Bluetooth® - Upper Arm Blood Pressure Monitor**

- Universal cuff: Also suitable for larger upper arm circumferences from 22-44 cm
- Risk indicator: Classification of measurements with coloured scale
- Ø morning and evening blood pressure for the last 7 days
- Wireless transmission of the measurement results





## **BR 60** Insect Bite Healer

- For treating insect bites and stings\*
- Also suitable if you are pregnant as there are no chemical additives - only by means of warmth
- Can help combat itching and swelling

\* Suitable for children over age 3  
(self-application from over age 12)





## **EM 49 Heat** Digital TENS/EMS Unit

- Transcutaneous electrical nerve stimulation (TENS)
- Electrical muscle stimulation (EMS)
- Relaxation and massage
- Electrode position display for easy programme selection and electrode placement
- 64 pre-programmed applications






## EM 59 Heat Digital TENS/EMS Device with Heat Function

- Transcutaneous electrical nerve stimulation (TENS)
- Electrical muscle stimulation (EMS)
- Relaxation and massage
- Electrode position display for easy programme selection and electrode placement
- Rechargeable battery





## **FT 15/1** Express Thermometer

- Contact-measurement technology
- Measuring time approx. 10 sec.
- Temperature alarm: Warning tone from 37.8°C
- Waterproof
- Can be disinfected
- Acoustic signal when measurement is complete





## FT 65 Multi-functional Thermometer

- 6-in-1 function: Ear, forehead and surface temperature, temperature alarm, date and time, 10 memory spaces
- Measurement in seconds
- Infrared measurement
- LED temperature alarm from 38°C





## **PO 30** Pulse Oximeter

- Determining oxygen saturation (SpO<sub>2</sub> ) and heart rate (pulse)
- Also suitable for sports at high altitudes
- Colour display with 4 available views
- Automatic switch-off

