



Always on hand for
insect bites and stings



OUTDOOR 
BR 10
With carabiner hook



BR 60
The family favourite



DAY & NIGHT 
BR 90
Incl. light

759.686_0624 Subject to errors and changes

Insect bite healers


Discover our entire product range
now at www.beurer.com

Beurer GmbH, P.O. Box 1427, 89004 Ulm,
Soeflinger Straße 218, 89077 Ulm, Germany,
Tel. +49 731/39 89-0. Fax. +49 731/39 89-139. export@beurer.de. www.beurer.com



*According to Medical Device Regulation (EU) 2017/745.
Observational study by the Prof. Dr. Kurscheid Institute to examine the effectiveness of the BR60 Insect
Bite Healer against pruritus and swelling after mosquito bites (MüPruritus) dated 16 February 2024 with
22 subjects. Further information can be found at www.beurer.com/br-study

beurer



Wave goodbye to itching

Our insect bite healers



Insect bite healer



The ideal insect bite healer for every scenario

A handy helper for outdoors –
now you can fight insect bites
and stings on the go



BR 10 Insect Bite Healer

Extra-small design with carabiner hook –
the perfect companion for everyday use,
when travelling and when exercising

- One programme (3 seconds)
- Product dimensions: 8.6 x 3.5 x 2.5 cm
- Medical device
- Incl. batteries



Insect bites are irritating

Mosquitoes, horseflies and wasps can often be unwelcome guests when swimming, barbecuing or exercising in the summer months. Insect bites and stings are difficult to avoid by natural means and can often have unpleasant consequences, such as a burning itching sensation and swelling. But our insect bite healers can help.

Principle behind the insect bite healer



Relief from itching and swelling

Use immediately after the insect bite or sting for the best results. The targeted heat can suppress itching and has been reported to encourage better healing. So easy to use, simply place the quick-heating ceramic hot plate on the affected area for as little as 3 seconds. Our insect bite healers do not use any chemicals, making them safe to use on children and in pregnancy.

The family favourite –
quick relief from itching
and pain



BR 60 Insect Bite Healer

Fits in any drawer or bag – for use at home or
on the move. Ideal for both adults and children

- 2 different programmes:
 - for initial use and persons with sensitive skin (3 seconds)
 - for regular use (6 seconds)
- Product dimensions: 9.5 x 3.5 x 2 cm
- Medical device
- Incl. batteries

A true all-rounder – hit the
target, even in the dark



BR 90 Insect Bite Healer

With optional light for precise treatment, even in the
dark. A real eye-catcher, with a stylish black design

- 2 different programmes:
 - for initial use and persons with sensitive skin (3 seconds)
 - for regular use (6 seconds)
- Product dimensions: 10.1 x 3.6 x 2 cm
- Medical device
- Incl. batteries