

# EM 50 Menstrual Relax


Many women struggle with period pain every month. Conditions such as endometriosis are also widespread and can cause unpleasant abdominal cramps. The EM 50 Menstrual Relax can offer pain relief – gently, quickly and without medication.

Comfortable to wear



Powerful battery

Optional heat function



## EM 50 Menstrual Relax

- Relief of menstrual pain, including in cases of endometriosis
- Combination of heat and electrical nerve stimulation (TENS)
- Heat and TENS can be applied separately
- Rechargeable battery
- Set includes: 1x Menstrual Relax device, 2x self-adhesive gel pads, 1x USB-C cable, instructions for use



With transcutaneous electrical nerve stimulation (TENS) to provide relief from period pain



15 TENS intensity levels




Heat function (43°C)

Discover our entire product range now at [www.beurer.com](http://www.beurer.com)

Beurer GmbH,  
P.O. Box 1427, 89004 Ulm, Soeflinger Strasse 218, 89077 Ulm, Germany.  
Tel. +49 731/39 89-0. Fax. +49 731/39 89-139.  
export@beurer.de [www.beurer.com](http://www.beurer.com)



759.684\_0722 Subject to errors and changes.



An end to menstrual pain!

Menstrual Relax




EM 50 Menstrual Relax  
TENS & Heat Pad



Pain relief through TENS and soothing heat



Practical for everyday use!  
The pad can easily be  
worn underneath clothing.



Feel better during  
your period

### Quick relief for abdominal pain

The EM 50 Menstrual Relax from Beurer combines TENS technology with soothing heat and offers women a whole new experience for managing abdominal pain. Thanks to the rechargeable battery and 15 intensity levels, you can decide for yourself where and when to use the pad!

### How does the EM 50 Menstrual Relax work?

The EM 50 Menstrual Relax targets the pain directly. It is placed on the lower abdomen with self-adhesive gel pads. Transcutaneous electrical nerve stimulation (TENS) suppresses the transmission of pain into the nerve fibres and relieves pain.



Ideal for use at home  
and on the move

### Soothing heat can be activated

The heat function of the EM 50 Menstrual Relax ensures that your body and muscles relax and that you feel completely comfortable, even on days of pain.

### Effective, without medication

By relieving pain with TENS, you can avoid taking painkillers, or at least reduce the amount that you need to take. After 20 minutes, the EM 50 Menstrual Relax will switch off automatically. Thanks to the rechargeable battery, the device is reusable.

### Flexible and discreet

The pad flexibly adapts to your body shape. It's also discreet, unobtrusive and comfortable to wear beneath your clothes when you're out and about.

