

Segmented display for precise view of individual regions of the body

PRECISE RESULTS AND SEGMENT ANALYSIS THANKS TO cross-measurement

Diagnostic bathroom scale for a new level of physical well-being

The BF 1000 Super Precision provides you with everything you need for comprehensive weight control. A precise full-body analysis can be generated using a special cross-measurement. The raised display on the aluminium stand and the sleek standing surface with ITO* coating allow for professional measurements at home!

*(ITO = Indium tin oxide)



"My experience shows that lifestyle changes through exercise and a changed diet are often not visible on conventional scales, which leads to people quitting. With bio-impedance scales, I see muscle development and fat loss, and remain motivated."

Prof. Dr. Thomas Kurscheid

BF 1000 Super Precision diagnostic bathroom scale

The BF 1000 Super Precision stands out with clever features, which make achievements as a result of exercise or weight loss very easy to measure.

- For professional weighing: raised dot matrix display on an aluminium stand with blue lighting for optimum readability
- Display of body weight, body fat percentage (incl. visceral fat), body water and muscle content, bone mass, calorie display, BMR and AMR
- Segmented display of individual body parts (arms, legs, torso) thanks to cross-measurement with 8 electrodes
- 10 user memories, each with 30 memory spaces + guest mode
- Connects to the free "beurer HealthManager" app
- Extra-sturdy: for 200 kg weight capacity with precise graduation of 50 g
- 5-year warranty




The app and software that are available for download may differ from the images. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. Android is a trademark of Google LLC.

Beurer GmbH,
P.O. Box 1427, 89004 Ulm, Soeflinger Strasse 218, 89077 Ulm, Germany
Tel. +49 731/39 89-0, Fax. +49 731/39 89-139,
export@beurer.de www.beurer.com www.beurer-healthguide.com



beurer



Stay attuned to all areas of your body – the professional scale for home use

BF 1000 SUPER PRECISION DIAGNOSTIC BATHROOM SCALE

For a detailed, full-body analysis



759.214_0719 Subject to errors and changes.



Why measure body fat?



Our weight is important for our health. The ratio of muscle and body fat percentage plays a decisive role, not just for athletes. Regular monitoring is primarily useful for the purpose of tracking achievements from exercise and weight loss, as well as avoiding becoming over or underweight and the conditions which result from this.

Monitoring success during fitness and weight training

Endurance athletes or weight lifters often want to have precise information about their body. In addition to just weight data, the status of the composition of muscle and fatty tissue is important for adjusting the workout to the athlete's current condition.

Monitoring success during a diet

Anyone who wishes to remain motivated during a diet will want to know how much fat their body has lost. Analysis of fat and muscle percentage can help to better adjust dietary plans on the path to a healthy level of physical well-being.

Is my body fat percentage healthy?

A slim figure is unfortunately not a guarantee for continuous good health. Visceral fat analysis (abdominal fat/organ fat) provides you with particularly important feedback, as even if someone appears to be healthy on the outside, fat cells can increasingly be accumulating in the centre of the body; this fat is linked to health risks.



Full-body analysis just like a professional

With the BF 1000 Super Precision you'll be informed of your weight values at all times. To perform the special cross-measurement, eight electrodes for the hands and feet have been integrated into the diagnostic bathroom scale. This produces results which are especially accurate.


In addition to body weight, the BF 1000 Super Precision can also determine your body fat, body water, muscle percentage and bone mass, and provide information about your calorie requirements.

Even visceral fat is displayed as well. Visceral fat surrounds the internal organs in the abdominal cavity. This abdominal/organ fat is suspected of contributing to conditions such as high blood pressure, vascular deposits or diabetes.

Weight 	Body fat 	Visceral fat 	Body water
Muscle percentage 	Bone mass 	BMI 	BMR/AMR (kcal)

More precise knowledge about your body than ever before

Additionally, the BF 1000 Super Precision offers a detailed segment analysis. It provides additional information about the fat and muscle percentage of individual body zones such as the left arm, right leg or torso, for example.




Display of fat and muscle percentage by body zones

Graphic evaluations in the modern dot matrix display allow your body values to be easily interpreted and consistently tracked via the app.



All values visible at a glance long-term

The BF 1000 Super Precision is compatible with the free "beurer HealthManager" app. Using Bluetooth®, all data can be transferred easily from the scale to a smartphone or tablet. This way you can keep an eye on your body values at all times, and easily reach your desired weight or training goals.



The "beurer HealthManager" app – a complete system for monitoring your health